



Implantable Lithium CFx

Features and Benefits

- Fully welded hermetically sealed design
- ISO 13485 certified
- High energy density
- Thermal shutdown separator
- Long shelf life
- Moderate rate capability

Applications

- Neuromodulation
- Monitoring
- Sleep apnea
- Incontinence
- Various
 - Rate and capacity dependent

Specifications

Part Number	95-1735
Weight	0.57g
Dimensions	See details on back
Open Circuit Voltage	3.4V
Cut-Off Voltage	2.0V
Capacity 37°C	65 mAh
Self-discharge 20°C	<1% Per year
Volume	0.273 cc
Operating Temperature	37°C (99°F)
Recommended Storage Conditions	20°C ± 5°C
Case Materials	Titanium
Case Polarity	Negative
Terminals	Titanium & Molybdenum
Transportation Safety	UN DOT ST/SG/AC.10/11/Revision
Battery Safety	IEC 60086-4
Options	Custom terminations Application specific testing

Contego P65 Battery

P65 37°C 4-Week Capacity Discharge
1 mA Pulse for 180 Seconds every 12 Hours
100 μ A Baseline, 70 mAh

