



Lithium CFx Battery

Features and Benefits

- Fully welded hermetically sealed design
- Case negative design
- Long shelf life
- Lightweight titanium case
- Proven chemistry

Applications

- Neuromodulators
- Cardiac
- Brain
- Stimulators

Lightweight, small cylindrical cell with excellent performance

Specifications

Part Number	95-1735
Weight	0.57 G
Dimensions	Height: 0.394" +/- 0.004" (10.0 mm) Diameter: 0.232" ± 0.002" (5.9 mm)
Open Circuit Voltage	3.4 V
Nominal Voltage	2.8 V
Cut-Off Voltage	2.0 V
Capacity 37 °C	65 mAh (28-Day rate)
Self-discharge 20°C	<1% Per year
Volume	0.273 cm ³
Energy Density	712 Wh/l
Specific Energy	347 Wh/kg
Operating Temperature	10°C – 55°C
Recommended Storage Conditions	20°C ± 5°C
Case Materials	Titanium
Case Polarity	Negative
Terminals	Titanium & Molybdenum
Transportation Safety	UN DOT ST/SG/AC.10/11/Revision 7 -15°C – 60°C

Typical performance based on historical results.

Contego P65 Battery

